Processing of Resterilizable Medical Devices

Zirc Dental Products, Inc.
3918 Hwy 55 SE
Buffalo, MN 55313 USA

Local: +1 763.682.6636

Fax: +1 763.682.6604

Products: Mouth Mirror, Mouth Mirror Heads and Mouth Mirror Handles.

EC REP Loser & Co. GmbH Benzstr. 1c, 51381 Leverkusen, Germany +49-2171-706670

NOTE: Please follow the instructions for use from the manufacturer of your cleaning and sterilization equipment. Hereinafter called “see manufacturer” information”.

1. **Restriction of Reprocessing:** Frequent reprocessing has marginal consequences for the Mouth Mirrors. Normally the end of durability of the product will be determined by visual degradation.
2. **Preparation Before Cleaning:** Place the Mouth Mirrors into your instrument holder (ex. Cassette).
3. **Cleaning:** There are three (3) cleaning methods that can be used to clean your Mouth Mirrors and Handles.
 - * Manual: Clean the Mouth Mirrors manually by washing with a non-corrosive, low sudsing neutral detergent.
 - * Ultrasonic: “see manufacturer information”
 - * Automated Washer: “see manufacturer information”After one of the above cleaning methods have taken place. Rinse Mouth Mirrors thoroughly.
4. **Drying:** Dry Mirrors thoroughly before sterilization. IMPORTANT: Failure to dry the mirrors thoroughly may cause marks to be burned into the glass surface or may lead to corrosion.
5. **Sterilization:** Steam sterilize for at least 5 minutes at 273° F/134° C or at least 20 minutes at 250° F/121° C. Do not heat above 277° F/136° C. To avoid corrosion the steam has to be free of ingredients. Recommended critical values of ingredients for table water and steam condensate are determined according to EN 285.

Please note: After sterilization in the Autoclave, please pay attention to the drying period INSIDE the Autoclave, do not take them out before drying period ends.

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