

Symbols Glossary

	Indicates the product manufacturer		Indicates the product distributor
	Indicates the product importer		Indicates the Authorized Representative in the European Union
	Indicates the Authorized Representative in Switzerland		Indicates the UK Responsible Person
	Country of Manufacture		Medical Device
	Product Catalogue Number		Product Lot Number
	Product Use-by Date		Product Date of Manufacture
	Identifies the unique device identifier		Keep Dry
	Keep away from sunlight		Indicates the temperature limits to which the product can be safely exposed.
	Do not use if package is damaged		Contains natural rubber or latex
	Non-sterile		Product is sterilizable in a steam sterilizer (autoclave)
	Product can be cleaned in ultrasonic bath		Product can be used with a washer-disinfector
	Consult Instructions for use or consult electronic instructions for use		Do not re-use
Rx ONLY	CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician		CE Registration Mark
	CE Registration Mark with notified body		Repackaging
	Repak Membership		Product does not contain latex