

Lingua-Fix®


Dry-Field Isolation

Made in the U.S.A.



Zirc

Because Time Is Everything



— Ø 6mm




- Placed Lingual Side
- Inserts directly into Saliva Ejector Valve
- Provides a dry work area
- Holds tongue in a steady position to protect from injuries
- Easy to place and remove
- Large suction surface avoids blockage and is more comfortable
- Reduces sedimentation of drilling by-products in pipes and tubing

10 Pack (50Z953S)

50 Pack (50Z953)

You might also like...
SE Cushion
(50Z970)



0413



Manufacturer
Zirc Dental Products, Inc.
Buffalo, MN 55313
USA | +1 763-682-6636
www.zirc.com

EC REP

Loser & Co. GmbH
Benzstr. 1C. D51381
Leverkusen, Germany
+49-2171-706670

20Z328-2 Rev M 0219