

STERI-SOAKER




DESCRIPTION

The Steri Soaker is an efficient, high-volume germicide tray designed for maximum use of chemicals to clean and disinfect/sterilize instruments. Designed for easy cleaning and is large enough to accommodate Zirc's Steri Containers, Cassettes, Bur Blocks and Bur Guard.

FEATURES

- Calendar on lid to remind staff when to change solution *(Use color code rings to mark dates)*
- Holds 1 gallon of liquid, fill line is marked in quarts and liters
- Safe-to-use chemicals; cold sterilants, cleaners (enzymes, detergents), disinfectants (phenols)
- To drain liquid from tray: lift from solution, move forward, set on grooves provided by base
- Hinged, fitted cover



Made in the U.S.A.

Zirc Dental Products, Inc.
3918 Hwy 55 SE, Buffalo, MN, USA 55313
+1 (763) 682-6636 • www.zirc.com



Losser & Co. GmbH Benzstr. 1C, 51381 Leverkusen, Germany +49-2171-706670



50Z805IFU REV A 0223



LEARN MORE ABOUT THIS PRODUCT